# CONTACTS

## AUSTRIA

Bauerfeind Ges.m.b.H.  
Hainburger Straße 33  
1030 Wien  
P +43 (0) 800 4430-130  
F +43 (0) 800 4430-131  
E info@bauerfeind.at

## BENELUX

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Dear customer,  
Thank you for choosing a Bauerfeind product.

We work to improve the medical effectiveness of our products every day – because your health is very important to us. Please closely observe the specifications in these instructions for use and the instructions given by the specialist.

If you have any questions, please contact your physician or medical retailer.

#### Intended purpose

SecuTec Genu Flex is a medical product. It is a stabilizing orthosis<sup>1</sup> that limits the range of motion in the case of complex knee injuries.

#### Indications

- Rupture of the anterior and / or posterior cruciate ligament (ACL / PCL)
- After ligament surgery / ligament reconstruction
- Severe and / or complex instability (traumatic, degenerative, e.g. "unhappy triad")
- Collateral ligament injuries
- To restrict the knee joint's range of motion (e.g. after meniscus refixation or meniscal implant)
- For non-surgical treatment or long-term rehabilitation
- Functional instabilities

#### Risks of using this product

- Only use the SecuTec Genu Flex according to your indication and the specific instructions given by a medical specialist<sup>2</sup>.
- Before using the product, the specialist or your physician must determine whether thrombosis prevention is necessary.
- Please consult a specialist or your physician if you want to use it together with other products. Do not modify the product yourself in any way; otherwise it may not help as expected or it may be harmful to your health. In these cases, any warranty or liability claims will be excluded.
- Avoid contact with substances containing grease or acids, ointments, lotions, extreme temperatures as well as chemicals with an oxidizing effect and strong acids.
- Minimize contact with abrasive substances (e.g. sand) in order to prevent stiff hinges.
- Any aids applied externally to the body may, if tightened excessively, lead to local pressure points or, in rare cases, constrict the underlying blood vessels or nerves.
- If you notice any changes or an increase in symptoms while wearing the product, stop any further use and contact your physician.
- Caution: When loosening the straps, or when removing the product or adjusting configured limitations of movement, there is not sufficient protection of the body part requiring treatment. In this case, particularly avoid placing any strain on this body part.
- If you become aware of any changes to the product, consult a medical specialist.
- If you want to use the SecuTec Genu Flex more than 24 months upon receipt, we recommend checking the plastic parts of the orthosis for damage prior to use because they can prematurely age due to sunlight, heat, cold and moisture, for example. Excessive use can limit the lifespan of the product.
- Should any problems with the fit of the orthosis arise, for instance if the shape of your leg changes (e.g. due to an increase or reduction in muscle), or should the product have a functional defect, please immediately contact your physician or attending specialist.
- Side effects involving the body as a whole have not been reported to date. Correct usage / fitting is assumed.
- Please adjust your behavior to the situation associated with your specific condition and any associated restriction of movement, and avoid anything that could delay recovery where possible. Proceed with particular caution.
- Furthermore, proceed with extreme caution, as your ability to stand and walk may be impacted.
- Ask your physician or the supervising specialist which specific restrictions are personally relevant for you and must be observed.
- If your perceptive faculties are restricted, if you take painkillers or blood-thinning medications or suffer from dementia, our products may only be worn when prescribed by a physician.

#### Contraindications

No clinically significant adverse reactions have been reported to date. If you have any of the following conditions, the product should only be used after consultation with your physician:

- Skin disorders / injuries to the relevant part of the body, particularly if inflammation is present. Likewise, any raised scars with swelling, redness, and excessive heat build-up
- Varicose veins (varicosis)
- Impaired sensation and circulatory disorders of the legs / feet (e.g. diabetes mellitus)

- Impaired lymphatic drainage – including unspecified soft tissue swelling located in distal direction of the fitted aid

#### Application instructions

##### Size determination:

The size is selected and fitted by a trained specialist as needed after measuring the patient. You can refer to the size chart on the packaging for the correct size.

##### Procedure for the first fitting of the SecuTec Genu Flex knee orthosis

The connecting element with the B logo defines the front side of the orthosis. It must be positioned on the front of the thigh. In addition, the upper end of the orthosis is marked by arrows on the frame. On Bauerfeind products, the arrows always point upwards.

##### Adjustment options for the hinge (only by a trained specialist)

- The hinge cover can be opened via a screw fastening. This cover fastens the movement-limiting wedges (extension – front / flexion – rear) and also protects the hinge mechanism against dirt.
- Opening the screw fastening (quarter turn until the arrow points to ) allows the hinge protection to be removed and the desired limiting wedges to be changed.
- The desired extension and flexion wedges are selected from the enclosed frame and are used  . The following limitations are possible:  
Extension: 0°, 10°, 20°, 30°, 45°  
Flexion: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- Then, click the hinge cover onto the hinge, check whether it is correctly positioned and firmly attached to the hinge, and close the screw fastening  . Make the adjustment to both hinges.

##### Adjustment of the frame (only by a trained specialist)

- To ensure an optimal fit, the hinges must be positioned parallel to the shin bone and the thigh. The middle joining elements are positioned centrally on the calf and the thigh.
- The circumference around the proximal and distal frame of the SecuTec Genu Flex can be adjusted by removing and adding the chain links.
- To do so, open the frame closure with a sharp object (e.g.: pen) .
- Separate the chain links from the cushion by pulling carefully on the cushion and thereby removing the plastic snap fasteners from the frame  a-b.
- Now you can loosen the connection between the chain links by bending against the curvature of the frame . To do so, grab the two chain links that you want to separate and bend them against the curvature of the frame by up to 90° to one another. Tip: First loosen the two outer guide edges by bending slightly, then the middle piece.
- As needed, you can now either remove the chain links or add new parts from the packaging included.
- Each chain link corresponds to a circumference change of 1 cm.
- To join together the chain links, first slide the two outer guide edges into the groove of the front chain link  a. After this, the center piece can be clicked in by bending slightly against the curvature of the frame  b.
- It is recommended to carry out alteration on the inside and outside of the orthosis symmetrically to the middle joining element.
- To reconnect the frame with the chain links, slide the frame closure into the frame mount .
- Finally, fasten the cushion back on the frame using the snap fasteners  a-b. Proceed symmetrically starting from the middle joining element.

##### Adjustment of the orthosis height (only by a trained specialist)

You can adjust SecuTec Genu Flex either distally or proximally in two lengths.

To do so, press a pointed object against the plastic button in the area of the frame mount  and slide the hinge leg upwards until it engages in the next position. Standard delivery is made in the shortest setting.

##### Crimping the hinge leg (only by a trained specialist)

If required, you can reshape the hinge leg using round crimping irons with protective caps. When doing so, the plastic parts can remain on the hinge legs.

##### Static fitting (only by a trained specialist)

- Perform the initial adjustment with the straps open.
- Bend the knee to be treated to a 30° angle and place the SecuTec Genu Flex on it so that the hinge pivot point of the orthosis hinge is aligned with the compromise axis of rotation as defined by Nietert .

##### Putting on the SecuTec Genu Flex

- Put the orthosis on against your skin. Wearing the support over pants can cause slipping and impair functionality.
- The leg is bent at a gentle angle (approximately 30°). The position of the kneecap is then determined through palpation.
- The center of the hinge (pivot point of the orthosis hinges) is positioned at approximately the same height as the center of the kneecap .
- Start by pulling the strap tight directly beneath the knee, then proceed to the straps closer to the joint above the knee. Then tighten the straps farther from the joint in the same order . When closing the straps, check whether the orthosis has rotated laterally around the leg and correct this if necessary.
- After all the straps have been closed, the orthosis must exert an even pressure when the muscles are tensed.
- Once the orthosis has been put on, check once more that it is in the correct position.

#### Part of the body these products are used for

In accordance with the indications (knee). Please see intended purpose.

#### Assembly and fitting instructions

The SecuTec Genu Flex is supplied as a pre-assembled orthosis in standard sizes. The first fitting, which includes the hinge adjustments along with the adjustment of the frame and orthosis height, crimping and static fitting, must be carried out by a trained specialist<sup>2</sup>.

#### Maintenance instructions

If handled and cared for correctly, the product is maintenance-free.

#### Technical specifications / parameters

The SecuTec Genu Flex is an orthosis for the knee joint. It consists of a frame with individual chain links, hinges, straps and cushions.

#### Notes on reuse

This product is intended to be fitted to one patient only.

#### Warranty

The statutory regulations of the country of purchase apply. Please first contact the retailer from whom you obtained the product directly in the event of a potential claim under the warranty. The product must be cleaned before submitting it in case of warranty claims. Warranty may be limited or excluded if the instructions on how to handle and care for the SecuTec Genu Flex have not been observed.

#### Warranty is excluded if:

- the product was not used according to the indication
- the instructions given by the specialist were not observed
- the product was modified arbitrarily

#### Notice of liability

For use of the SecuTec Genu Flex, the following applies:

#### Taking off the SecuTec Genu Flex

- When taking off the orthosis, it is recommended to adopt a seated position.
- Open the straps marked in the image  in reverse order and slide the orthosis down on the leg.  
TIP: If you have been strictly instructed not to bend your leg, ask someone around you to help pull the orthosis over your foot.

#### Cleaning instructions

Never expose the orthosis to direct heat (e.g. heater, sunlight, etc.)!

- The aluminum frames of the orthosis are coated and can be cleaned using pH-neutral soap like the plastic parts.
- We recommend washing the cushions and straps by hand at 30 °C.
- In case of contact with salt water or chlorinated water, thoroughly rinse the orthosis out with plain water afterwards.
- If you notice that a hinge is stiff after contact with coarse dirt (e.g. with sand), contact your specialist to have the hinge cleaned. This product has been tested by our own quality control system. However, should you have any complaints, please contact your local medical supply retailer.

Do not diagnose yourself or decide on which medication to take unless you are a medical specialist. Before first using our product, please seek the advice of this medical specialist, as this is the only way to assess the effect of our SecuTec Genu Flex on your body and determine the potential risks of using this product due to your personal constitution. Follow the advice of the medical specialist as well as all instructions in this brochure or its online version, including any excerpts (as well as texts, images, graphics etc.). If you still have doubts after your consultation with the specialist, please contact your physician or retailer, or contact us directly.

#### Duty to report

Due to regional legal regulations, you are required to immediately report any serious incident involving the use of this medical device to both the manufacturer and the responsible authority. Our contact details can be found on the back of this brochure.

#### Disposal

Leave orthosis: Once you have finished using the product, please dispose of it according to local specifications or return the product back to the medical retailer who supplied you with the orthosis.

#### Material content

Polyketone (PK), Aluminum (AL), Polyamide (PA), Polyurethane (PUR), steel, polished, anodized, Polyester (PES), Polyoxymethylene (POM), Thermoplastic polyurethane (TPU), Polypropylene (PP), Elastane (EL), stainless steel, Glass fiber (GF), Polyethylene terephthalate (PET), brass, nickel-free, adhesive gum resin

 – Medical device

 – DataMatrix Unique Device Identifier UDI

<sup>1</sup> Orthosis = orthopedic aid used to stabilize, relieve, immobilize, control, or correct the limbs or the torso

<sup>2</sup> A specialist is any person who is authorized, according to their state's regulations, to fit and provide instructions on the use of supports and orthoses.

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